



ResMed

Actual Size:
148 mm (Width)
x 105 mm (Height)



Ruler is for scaling purposes.


How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

AirFit™ N30i

Nasal cradle mask







ResMed

AirFit™ N30i
Nasal mask

AirFit™ P30i
Nasal pillows mask

Frame Fitting Instructions



 Small frame	 Standard frame
0-150mm	>150mm
0-6"	>6"